



#### INTRODUCTION



In the business of Manufacturing & Global Distribution of Speciality chemicals

R&D to create & Produce niche products in India

Corporate office in Mumbai.
Manufacturing at Tarapur ISO
13485:2003 certified


### THE PRODUCT & FEATURES



Simply place in

enclosed region. Body

Parts rich in blood



Designed for troops HAA/SHAA sectors to overcome extreme cold related afflictions




3 years of extensive R& D and certified by NABL accredited laboratory from Bangalore

like Chill Blains, Frost bite and Hypothermia induced by weather conditions

#### THE PRODUCT & FEATURES


No fumes.






## THE PRODUCT & FEATURES







## TEMPERATURE CURVE OF W-ARMY










Increases blood circulation and recovery of normal physiological function

Enhances efficiency, combat readiness and psychological morale







"Fools" Glomus Bodies and counters body aberration



Reduces the burden on the existing medical infrastructure



### W-ARMY IS THE RIGHT FIT FOR TROOPS



Undergoing acclimatization

LR-CI Patrolling in low temperature & snow bound areas Posts above 9000 Ft

Med-evac Cases



#### W-ARMY IS THE RIGHT FIT FOR TROOPS



Search and Rescue operations

Air crew and passengers flying in non pressurized aircraft and heptors Reconnaissance expeditions

## OUR PRODUCT RANGE



Regular



Foot warmer



Mini hand warmer



Back-Shoulder-Neckstomach




Toe Pair




Sticking type



## PRODUCT USAGE



Regular



Mini



Toe



Foot Sole



# SPECIFICATIONS

2711				SAL VIII GET P			DATE AND 4	ACCUPATION OF THE PARTY OF THE	
TYPE/SHAPE	MODEL NO.	PACKING Pouch per Pack	SIZE MM (+/- 5mm) (LENGTH* WIDTH)	SIZE TOLERANCE	INNER PACK WEIGHT (gm)	BEAKER TEST for MIXTURE	TIME TO REACH 45 °C (MINUTES)	MAX TEMP°C	DURATIO N 45°CTO 45°C HR
REGULAR	WAEE01		120 * 80	+/- 5MM	60 + / - 5	MIN 80 °C IN 40 MIN	15 TO 20	64 - 65	MIN 8 HRS
MINI	WAEE02		90 * 70	+/- 5MM	30 + / - 3	MIN 80 °C IN 40 MIN	15 TO 20	55 - 56	MIN 8 HRS
TOE PAIR	WAEE03	2	90 * 70	+/- 5MM	15+/-2	MIN 80 °C IN 40 MIN	15 TO 20	48- 52	MIN 5 HRS
FOOT PAIR	WAEE04	2	250 * 75	+/-10MM	30+/ - 2	MIN 80 °C IN 40 MIN	15 TO 20	45 - 48	MIN 5 HRS
SHOULDER	WAEE05	1	250 * 90	+/-10MM	45 +/-3	MIN 80 °C IN 40 MIN	15 TO 20	50 - 55	MIN 5 HRS
REGULAR STICK	WAES01	1	120 * 80	+/- 5MM	40 + / -3	MIN 80 °C IN 40 MIN	15 TO 20	55 - 60	MIN 8 HRS
MINI - STICK	WAES02	7	90 * 70	+/- 5MM	20 + / - 2	MIN 80 °C IN 40 MIN	15 TO 20	50 - 55	MIN 8 HRS





#### WAY FORWARD

Scope for expanding the benefits of W-Army to a larger user base



Significant user appreciation – Army Recommendation

Ready to conduct trials and demos at any location if required Imperative part of standard gear for Trps in HAA/SHAA. Help increase operational efficiencies for peace time & Battle/Combat missions...

